

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/GREY/CHILD (S)

Réf. 330021



VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/GREY/CHILD (S)



Réf. 330021

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446006822
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	15,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,138
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference 330021

All characteristics are given for information only and may be changed