

VAQUEZ-LAUBRY NANO/2 TUBES/ADULT BLACK NYLON CUFF (M)

Réf. 330210



VAQUEZ-LAUBRY NANO/2 TUBES/ADULT BLACK NYLON CUFF (M)



Réf. 330210

Made of black Oxford-grade nylon, with inner bladder and dual tubing (10 cm & 60 cm), with split leather manometer holder, and Velcro fastening.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446028596
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,180
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference 330210

All characteristics are given for information only and may be changed