

LIAN CLASSIC COTTON CUFF WITH BLADDER/2 TUBES/GREY/LARGE ADULT (L)



Réf. 332001



LIAN CLASSIC COTTON CUFF WITH BLADDER/2 TUBES/GREY/LARGE ADULT (L)



Réf. **332001**

Cotton cuff with highly resistant rubber bladder.
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
Washable at 30-40°C.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446028466
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	7,1
Length (cm) :		Length packaging (cm) :	14,3
Height (cm) :		Height packaging (cm) :	0,4
Weight (kg) :		Weight packaging (kg) :	0,310
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:23 - This article is available on our website : www.gsh-med.fr with the following reference 332001

All characteristics are given for information only and may be changed