

## LIAN CLASSIC COTTON CUFF WITH BLADDER/2TUBES/CARTOON DESIGN/CHILD (S)

Réf. 332029



# LIAN CLASSIC COTTON CUFF WITH BLADDER/2TUBES/CARTOON DESIGN/CHILD (S)



Réf. 332029

Cotton cuff with highly resistant rubber bladder.  
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
Washable at 30-40°C.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446008796
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,094
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 332029

All characteristics are given for information only and may be changed