


## LIAN CLASSIC COTTON CUFF W. BLADDER/2TUBES/CARTOON DSGN/NEWBORN (XXS)

Réf. 332049



## LIAN CLASSIC COTTON CUFF W. BLADDER/2TUBES/CARTOON DSGN/NEWBORN (XXS)



Réf. 332049

Cotton cuff with highly resistant rubber bladder.  
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
Washable at 30-40°C.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446008772
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	13,5
Height (cm) :		Height packaging (cm) :	2,5
Weight (kg) :		Weight packaging (kg) :	0,034
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:23 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 332049

All characteristics are given for information only and may be changed