

VAQUEZ-LAUBRY CLINIC CUFF/LARGE ADULT (L)

Réf. 338000



VAQUEZ-LAUBRY CLINIC CUFF/LARGE ADULT (L)



Réf. 338000

Made of yellow Polymex®, with built-in obturating system to allow decontamination by immersion, and Velcro fastening. Legal markings engraved and non-removable.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446006594
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	18,0
Length (cm) :		Length packaging (cm) :	24,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,126
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:35 - This article is available on our website : www.gsh-med.fr with the following reference 338000

All characteristics are given for information only and may be changed