

## STETHOSCOPE MAGISTER ADULTE DOUBLE PAVILLON GRIS ALIZE



Réf. 507204



CE



# STETHOSCOPE MAGISTER ADULTE DOUBLE PAVILLON GRIS ALIZE



**Réf. 507204**

Pearly tube model. For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446004026
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,398
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/12/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507204

All characteristics are given for information only and may be changed