

SPENGLER PAEDIATRIC MAGISTER STETHOSCOPE, DUAL, BLACK

Réf. 507210



SPENGLER PAEDIATRIC MAGISTER STETHOSCOPE, DUAL, BLACK



Réf. 507210

Stethoscope, binaural, complete: A mechanical listening device designed for listening to sounds from the heart and lungs from paediatric patients.

For paediatric auscultation.

- Double stainless steel chestpiece: adult and a bell chestpiece
 - Dual frequency membrane
 - Light pressure to hear low frequencies (vascular and cardiac)
 - Firm pressure to hear high frequencies (pulmonary)
 - Light stainless steel binaural, designed for the auricular anatomy
 - Comfortable ear tips and non-chill ring
- Black colour.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009120
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	16,5
Length (cm) :		Length packaging (cm) :	32,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,375
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/12/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference 507210

All characteristics are given for information only and may be changed