

MAGISTER STETHOSCOPE FOR CHILDREN / SKY BLUE

Réf. 507214



MAGISTER STETHOSCOPE FOR CHILDREN / SKY BLUE



Réf. 507214

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009137
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,375
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



01/12/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference 507214

All characteristics are given for information only and may be changed