

VAQUEZ-LAUBRY CLINIC / IN POLYBAG / LARGE ADULT (L)

Réf. 514200



VAQUEZ-LAUBRY CLINIC / IN POLYBAG / LARGE ADULT (L)



Réf. 514200

Cuff-mounted designed to effectively prevent nosocomial infections. The manometer attached on the cuff can be easily removed for liquid decontamination of the cuff in most disinfection solutions on the market. Large, shock-resistant ABS dial and large inflation bulb. 3 cuff sizes available. Comes in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446006617
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,468
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:28 - This article is available on our website : www.gsh-med.fr with the following reference 514200

All characteristics are given for information only and may be changed