

VAQUEZ-LAUBRY CLINIC / IN POLYBAG / ADULT (M)

Réf. 514201



VAQUEZ-LAUBRY CLINIC / IN POLYBAG / ADULT (M)



Réf. 514201

Cuff-mounted designed to effectively prevent nosocomial infections. The manometer attached on the cuff can be easily removed for liquid decontamination of the cuff in most disinfection solutions on the market. Large, shock-resistant ABS dial and large inflation bulb. 3 cuff sizes available. Comes in a nylon pouch.

| FEATURES | | | |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device | Yes | Classe Im (fonction de mesurage) | |
| Marquage | CE0459 | Barcode | 3700446008604 |
| Guarantee (month): | 24 | Customs Code | 90189010 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm) : | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,434 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -20 | Temperature max | 70 |
| Humidity min | 20 | Humidity max | 85 |



31/03/2022 02:28 - This article is available on our website : www.gsh-med.fr with the following reference 514201

All characteristics are given for information only and may be changed