

ADULT CUFF (M) FOR MOBI

Réf. 513712



ADULT CUFF (M) FOR MOBI



Réf. 513712

| FEATURES | | | |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device | Yes | Classe Im (fonction de mesurage) | |
| Marquage | | Barcode | 3700446028015 |
| Guarantee (month): | 12 | Customs Code | 90189010 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,276 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



26/04/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference 513712

All characteristics are given for information only and may be changed