

## VAQUEZ-LAUBRY CLINIC / MULTI CUFF / (L, M, S)

Réf. 514300



## VAQUEZ-LAUBRY CLINIC / MULTI CUFF / (L, M, S)



**Réf. 514300**

Cuff-mounted designed to effectively prevent nosocomial infections. The manometer attached on the cuff can be easily removed for liquid decontamination of the cuff in most disinfection solutions on the market. Large, shock-resistant ABS dial and large inflation bulb. 3 cuff sizes available. Comes in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024765
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,846
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:28 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 514300

All characteristics are given for information only and may be changed