

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/ADULT (M)/GREY

Réf. 518011



# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/ADULT (M)/GREY



Réf. 518011

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446007287
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	8,0
Weight (kg) :		Weight packaging (kg) :	0,750
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:13 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 518011

All characteristics are given for information only and may be changed