

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE

Réf. 518002



# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



**Réf. 518002**

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446007294
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,831
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:28 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 518002

All characteristics are given for information only and may be changed