

CUFF/2 TUBES (60+10CM)/NYLON/WITH BULB/GREY/LARGE ADULT (L)

Réf. 335101



CUFF/2 TUBES (60+10CM)/NYLON/WITH BULB/GREY/LARGE ADULT (L)



Réf. 335101

Cuffs with dual tube bladder & coupling (60 cm tube with Classic® bulb and valve, and 10 cm with male screw-on metal fitting)

FEATURES			
Medical device	Yes	Classe I	
Marquage		Accuracy measurement	Classe I
Guarantee (month):	6	Barcode	3700446028459
Sterile	No	Customs Code	9018901000
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,310
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 335101

All characteristics are given for information only and may be changed