

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)

Réf. DS99306



CE
0483



ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Réf. **DS99306**

Adult cuff (26 to 35 cm)
Washable nylon cuff without pocket with Velcro fastening

| FEATURES | | | |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device | Yes | Classe Im (fonction de mesurage) | |
| Marquage | CE0483 | Barcode | 3661040075044 |
| Guarantee (month): | 12 | Customs Code | 9018901000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 8,8 | Width packaging (cm): | |
| Length (cm) : | 16,0 | Length packaging (cm) : | |
| Height (cm) : | 4,0 | Height packaging (cm) : | |
| Weight (kg) : | 0,093 | Weight packaging (kg) : | 0,093 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -20 | Temperature max | 70 |
| Humidity min | 10 | Humidity max | 85 |



30/03/2022 14:59 - This article is available on our website : www.gsh-med.fr with the following reference DS99306

All characteristics are given for information only and may be changed