

EXTRA LARGE CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(32 TO 42 CM)

Réf. DS99606



EXTRA LARGE CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(32 TO 42 CM)



Réf. **DS99606**

Extra large cuff (32 to 42 cm)

Washable nylon cuff without pocket with Velcro fastening

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0483	Barcode	3661040075310
Guarantee (month):	12	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	
Length (cm) :	27,2	Length packaging (cm) :	
Height (cm) :	6,8	Height packaging (cm) :	
Weight (kg) :	0,122	Weight packaging (kg) :	0,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	10	Humidity max	85



30/03/2022 14:59 - This article is available on our website : www.gsh-med.fr with the following reference DS99606

All characteristics are given for information only and may be changed