

## CHILD CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(18 TO 26 CM)



Réf. DS99406



CE  
0483



# CHILD CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(18 TO 26 CM)



Réf. **DS99406**

Child cuff (18 to 26 cm)  
Washable nylon cuff without pocket with Velcro fastening

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0483	Barcode	3661040075211
Guarantee (month):	12	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,0	Width packaging (cm):	
Length (cm) :	12,2	Length packaging (cm) :	
Height (cm) :	4,6	Height packaging (cm) :	
Weight (kg) :	0,058	Weight packaging (kg) :	0,058
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	10	Humidity max	85



30/03/2022 14:59 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DS99406

All characteristics are given for information only and may be changed