

## SOFT EARPIECE FOR IDEAL+ (PAIR), RED



Réf. DS8009\_R



## SOFT EARPIECE FOR IDEAL+ (PAIR), RED



Réf. DS8009\_R

FEATURES			
Medical device	No		
Marquage		Barcode	3661040073972
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,0	Width packaging (cm):	
Length (cm) :	1,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,006	Weight packaging (kg) :	0,004
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:16 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DS8009\_R