

COOL BAG FOR DIABETIQUE PM DIA, BLUE

 ELITE BAGS

Réf. MEDIA02_B



COOL BAG FOR DIABETIQUE PM DIA, BLUE



Réf. MEDIA02_B

FEATURES			
Medical device	No		
Marquage	Barcode	8436564540690	
Guarantee (month):	24	Customs Code	4202.12.99.90
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,0	Width packaging (cm):	12,0
Length (cm) :	19,0	Length packaging (cm) :	20,0
Height (cm) :	5,0	Height packaging (cm) :	4,5
Weight (kg) :	0,130	Weight packaging (kg) :	0,214
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



24/11/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference MEDIA02_B