

## COOL BAG FOR DIABETIQUE PM DIA, GREY

 ELITE BAGS

Réf. MEDIAB03

GARANTIE  
**2**  
ANS  
YEARS  
GUARANTEED

CE



## COOL BAG FOR DIABETIQUE PM DIA, GREY



Réf. **MEDIAB03**

FEATURES			
Medical device	No		
Marquage	CE	Barcode	8436564542199
Guarantee (month):	24	Customs Code	4202.12.99.90
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,0	Width packaging (cm):	19,0
Length (cm) :	17,0	Length packaging (cm) :	10,5
Height (cm) :	6,0	Height packaging (cm) :	6,0
Weight (kg) :	0,012	Weight packaging (kg) :	0,248
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	30	Humidity max	85



24/11/2022 21:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference MEDIAB03