

NANCY_T00

Réf. WSMED07_00



NANCY_T00



Réf. WSMED07_00

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040197883
Guarantee (month):	12	Customs Code	62114310030
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	30,0	Width packaging (cm):	30,0
Length (cm) :	38,0	Length packaging (cm) :	38,0
Height (cm) :	1,5	Height packaging (cm) :	1,5
Weight (kg) :	0,300	Weight packaging (kg) :	0,300
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:04 - This article is available on our website : www.gsh-med.fr with the following reference WSMED07_00