

LIAN NM / GREY NYLON CUFF / LARGE ADULT (L)

Réf. 520620



CE
0459



LIAN NM / GREY NYLON CUFF / LARGE ADULT (L)



Réf. 520620

Ambidextrous - central decompression valve. Non-deformable, chromium-plated metal casing with foam manometer cover. Large inflation bulb, 6 cuff sizes available with metal quick coupling. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446024796
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,736
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	85	Humidity max	85



26/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference 520620

All characteristics are given for information only and may be changed