

## FEMALE TUNIC, IRIS LYOCELL , S3, HOLTEX



Réf. WMIRI07\_93



## FEMALE TUNIC, IRIS LYOCELL , S3, HOLTEX



Réf. WMIRI07\_93

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661041009437
Guarantee (month):	0	Customs Code	6211433100
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference WMIRI07\_93