

Réf. MEEMER04_B

GARANTIE
2
ANS
YEARS
GARANTEE



Réf. MEEMER04_B

Bag designed to carry health equipment linked to oxygen therapy.

- Internal Subdivision in two heights:
 - top removable tray with adjustable dividers,
 - lower part with three separators, to store the oxygen bottle, and a mesh pocket in the flap.
- It also has three external pockets, including a front one with an internal subdivision with additional pockets and elastic bands.
- High resistance zips
- Reflective strips
- Reinforced handles and straps
- External card holder
- Transparent top pocket
- Dimensions 63 x 22.5 x 24 cm; capacity 28 l.
- Material 600D Polyester
- Weight 1.89 kg

FEATURES			
Medical device	No		
Marquage		Barcode	3661040001180
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	24,0	Width packaging (cm):	23,0
Length (cm) :	63,0	Length packaging (cm) :	63,0
Height (cm) :	22,5	Height packaging (cm) :	24,0
Weight (kg) :	1,890	Weight packaging (kg) :	1,890
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER04_B