

## MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 40, GREEN

Réf. WSAM\_40V

medi Plogs



CE

DISPOSITIF MEDICAL  
MEDICAL DEVICE

CLASSE  
I



# MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 40, GREEN



Réf. WSAM\_40V

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040100036
Guarantee (month):	0	Customs Code	64029993
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	
Length (cm) :	31,0	Length packaging (cm) :	
Height (cm) :	11,0	Height packaging (cm) :	
Weight (kg) :	1,000	Weight packaging (kg) :	1,160
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference WSAM\_40V