

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 42, NAVY BLUE

Réf. WSAM_42BF

medi Plogs



CE

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
I



MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 42, NAVY BLUE



Réf. WSAM_42BF

Clogs suitable for modern hospitals, operating theatres and intensive care units.
Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040100128
Guarantee (month):	0	Customs Code	64029993
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	22,0
Length (cm) :	31,0	Length packaging (cm) :	31,0
Height (cm) :	11,0	Height packaging (cm) :	12,0
Weight (kg) :	1,000	Weight packaging (kg) :	0,850
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:36 - This article is available on our website : www.gsh-med.fr with the following reference WSAM_42BF