

MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 45, FUCHSIA

Réf. WSAM_45L

medi Plogs



CE

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
I



MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 45, FUCHSIA



Réf. WSAM_45L

Clogs suitable for modern hospitals, operating theatres and intensive care units.
Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

E.g : Ref. WSAM45L for a pair of white size 45Mediplogs clogs

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040100319
Guarantee (month):	0	Customs Code	64029993
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	
Length (cm) :	31,0	Length packaging (cm) :	
Height (cm) :	11,0	Height packaging (cm) :	
Weight (kg) :	1,000	Weight packaging (kg) :	1,160
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:37 - This article is available on our website : www.gsh-med.fr with the following reference WSAM_45L

All characteristics are given for information only and may be changed