

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 47, NAVY BLUE

Réf. WSAM_47BF

medi Plogs



CE



MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 47, NAVY BLUE



Réf. WSAM_47BF

Clogs suitable for modern hospitals, operating theatres and intensive care units.
Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | CE | Barcode | 3661040100425 |
| Guarantee (month): | 0 | Customs Code | 64029993 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 22,0 | Width packaging (cm): | 22,0 |
| Length (cm) : | 31,0 | Length packaging (cm) : | 31,0 |
| Height (cm) : | 11,0 | Height packaging (cm) : | 12,0 |
| Weight (kg) : | 1,000 | Weight packaging (kg) : | 1,150 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



14/10/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference WSAM_47BF