

MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 47, WHITE

Réf. WSAM_47W

mediPlogs



CE



MEDIPLUGS OPERATING ROOM SHOES, PERFORATED, 47, WHITE



Réf. WSAM_47W

- Clogs suitable for modern hospitals, operating theatres and intensive care units.
Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
 - The back of the heel is sloped and rounded to make walking easier
 - Ventilation system allows air to circulate
 - Toe prints to enhance hold and spurs to aid the circulation
 - 100% CFC-free polyurethane
 - Washable by hand or machine at 60°C
 - Sizes: 35 - 47
- 5 Colors: white, fuschia,navy blue, green, light blue

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040100463
Guarantee (month):	0	Customs Code	64029993
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	
Length (cm) :	31,0	Length packaging (cm) :	
Height (cm) :	11,0	Height packaging (cm) :	
Weight (kg) :	1,000	Weight packaging (kg) :	1,160
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference WSAM_47W

All characteristics are given for information only and may be changed