

AJUSTABLE STOOL WITH NYLON WHEELS5 DOUBLE CASTORS BLACK



Réf. OS007_N



AJUSTABLE STOOL WITH NYLON WHEELS5 DOUBLE CASTORS BLACK



Réf. OS007_N

- Adjustable Stool nylon wheels 5 feet sitting 3.5cm
- 1 5-arm structure ABS diameter 58 cm
 - 5 swivel castors clip
 - 1 mini piston height adjustable gas 53 cm / max 70 cm
 - 1 seat foam 30 kg / M3 diameter 43 cm with screwed steel plate
 - Washable plastic coating retardant M1 classification
 - Load: 150 kg

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040162645
Guarantee (month):	24	Customs Code	9402900000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	17,5
Length (cm) :		Length packaging (cm) :	41,0
Height (cm) :		Height packaging (cm) :	42,0
Weight (kg) :	5,200	Weight packaging (kg) :	6,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:40 - This article is available on our website : www.gsh-med.fr with the following reference OS007_N

All characteristics are given for information only and may be changed