

## VAQUEZ LAUBRY NANO VELCRO CUFF-BLACK NYLON ADULT.LARGE SPENGLER

Réf. 512200



# VAQUEZ LAUBRY NANO VELCRO CUFF-BLACK NYLON ADULT.LARGE SPENGLER



Réf. 512200

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446007263
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	8,0
Weight (kg) :		Weight packaging (kg) :	0,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 01:22 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 512200

All characteristics are given for information only and may be changed