

## LITESTICK PENLIGHT - BLACK

Réf. 570620



## LITESTICK PENLIGHT - BLACK



Réf. **570620**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009830
Guarantee (month):	24	Customs Code	90189080907
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,046
SPECIFIC STORAGE CONDITION			
Temperature min	10	Temperature max	40
Humidity min	(N/A)	Humidity max	(N/A)



31/08/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 570620