

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK

Réf. 507300



CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Réf. 507300

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass.
Elegant and streamlined for easy handling.
High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips.
Ideal for monitoring blood pressure and general auscultation.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009151
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	16,5
Length (cm) :		Length packaging (cm) :	32,5
Height (cm) :		Height packaging (cm) :	0,5
Weight (kg) :		Weight packaging (kg) :	0,393
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : www.gsh-med.fr with the following reference 507300

All characteristics are given for information only and may be changed