

CHROME-PLATED LAUBRY STETHOSCOPE / NAVY BLUE

Réf. 507302



CHROME-PLATED LAUBRY STETHOSCOPE / NAVY BLUE



Réf. 507302

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009168
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,393
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : www.gsh-med.fr with the following reference 507302

All characteristics are given for information only and may be changed