

## CHROME-PLATED LAUBRY STETHOSCOPE / RED

Réf. 507308



# CHROME-PLATED LAUBRY STETHOSCOPE / RED



**Réf. 507308**

The famous cardiophone. Single head chestpiece ( in Zamak - 107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446026523
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,393
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507308

All characteristics are given for information only and may be changed