

## DUAL PULSE STETHOSCOPE / NAVY BLUE

Réf. 507422



## DUAL PULSE STETHOSCOPE / NAVY BLUE



Réf. **507422**

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009212
Guarantee (month):	36	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	24,5
Length (cm) :		Length packaging (cm) :	11,5
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,237
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:45 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507422

All characteristics are given for information only and may be changed