

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)

Réf. 512220



CE
0459



VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Réf. 512220

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0459	Barcode	3700446007256
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,582
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:27 - This article is available on our website : www.gsh-med.fr with the following reference 512220

All characteristics are given for information only and may be changed