

## LIAN NM / NAVY BLUE COTTON CUFF / ADULT (M)

Réf. 520627



CE  
0459



# LIAN NM / NAVY BLUE COTTON CUFF / ADULT (M)



Réf. 520627

Ambidextrous - central decompression valve. Non-deformable, chromium-plated metal casing with foam manometer cover. Large inflation bulb, 6 cuff sizes available with metal quick coupling. Supplied in a nylon pouch.

| FEATURES                   |           |                                  |               |
|----------------------------|-----------|----------------------------------|---------------|
| Medical device             | Yes       | Classe Im (fonction de mesurage) |               |
| Marquage                   | CE0459    | Barcode                          | 3700446008864 |
| Guarantee (month):         | 36        | Customs Code                     | 90189010      |
| Sterile                    | No        |                                  |               |
| INFORMATIONS               | Main unit | Packaged unit                    |               |
| Width (cm) :               |           | Width packaging (cm):            |               |
| Length (cm) :              |           | Length packaging (cm) :          |               |
| Height (cm) :              |           | Height packaging (cm) :          |               |
| Weight (kg) :              |           | Weight packaging (kg) :          | 0,691         |
| SPECIFIC STORAGE CONDITION |           |                                  |               |
| Temperature min            | -20       | Temperature max                  | 70            |
| Humidity min               | 85        | Humidity max                     | 85            |



26/04/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 520627

All characteristics are given for information only and may be changed