

## LARGE ADULT CUFF (L) FOR MAXI+3 SPENGLER

Réf. 522712



# LARGE ADULT CUFF (L) FOR MAXI+3 SPENGLER



Réf. **522712**

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage		Barcode	3700446028183
Guarantee (month):	60	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,118
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:20 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 522712

All characteristics are given for information only and may be changed