

CHILD CUFF (S) FOR MAXI+3 SPENGLER



Réf. 522713



CHILD CUFF (S) FOR MAXI+3 SPENGLER



Réf. 522713

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage		Barcode	3700446028190
Guarantee (month):	60	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,071
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:20 - This article is available on our website : www.gsh-med.fr with the following reference 522713

All characteristics are given for information only and may be changed