


LED FIBER OPTIC OTOSCOPE - SMARTLED 5500

Réf. 570517



LED FIBER OPTIC OTOSCOPE - SMARTLED 5500



Réf. **570517**

LED lighting with Optical fibre
 The Spengler otoscope SMARTLED 5500® is using LED lighting with optical fibre allowing a direct, homogeneous illumination and without reflections of the auricular canal and the eardrum. Luminance 85K LUX. Bimaterial head in ABS and Stainless steel, screwed head for a safe fixation. Ergonomic ABS grip. End clamp with built-in switch button. Optical fibre for optimal transmission of the light beam. x 3 magnification. Supplied in a case with two LR6 AA batteries, 8 speculums (two of 2.5 mm and four of 4 mm).

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446026332
Guarantee (month):	36	Customs Code	90189042
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	11,5
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,385
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/05/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference 570517

All characteristics are given for information only and may be changed