

## LITESTICK PENLIGHT - SILVER

Réf. 570624



## LITESTICK PENLIGHT - SILVER



Réf. 570624

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009854
Guarantee (month):	24	Customs Code	901890
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,046
SPECIFIC STORAGE CONDITION			
Temperature min	10	Temperature max	40
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:54 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 570624

All characteristics are given for information only and may be changed