




## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 BLUEBERRY



Réf. 570517B



## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 BLUEBERRY



Réf. 570517B

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446041717
Guarantee (month):	36	Customs Code	90189042
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	20,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,385
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	10	Humidity max	90



31/05/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 570517B

All characteristics are given for information only and may be changed