




LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 CORAL



Réf. 570517C



LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 CORAL



Réf. **570517C**

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | CE | Barcode | 3700446041700 |
| Guarantee (month): | 36 | Customs Code | 90189042 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | 12,0 |
| Length (cm) : | | Length packaging (cm) : | 20,0 |
| Height (cm) : | | Height packaging (cm) : | 4,0 |
| Weight (kg) : | | Weight packaging (kg) : | 0,385 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | -20 | Temperature max | 70 |
| Humidity min | 10 | Humidity max | 90 |



31/05/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference 570517C

All characteristics are given for information only and may be changed