

## DORO MAGNA 4000, 60DB



Réf. HDMAGN01N



## DORO MAGNA 4000, 60DB



Réf. HDMAGN01N

FEATURES			
Medical device	No		
Marquage	Barcode	7322460063795	
Guarantee (month):	24	Customs Code	85171800
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	10,0
Length (cm) :	22,5	Length packaging (cm) :	23,0
Height (cm) :	9,8	Height packaging (cm) :	18,0
Weight (kg) :	1,150	Weight packaging (kg) :	1,208
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference HDMAGN01N