

DORO PHONE CORDLESS DECT PHONEEASY 100W DUO, WHITE



Réf. HDPHON13W



DORO PHONE CORDLESS DECT PHONEEASY 100W DUO, WHITE



Réf. **HDPHON13W**

FEATURES			
Medical device	No		
Marquage		Barcode	7322460055493
Guarantee (month):	24	Customs Code	8517110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,0	Width packaging (cm):	10,0
Length (cm) :	19,0	Length packaging (cm) :	19,0
Height (cm) :	18,0	Height packaging (cm) :	18,0
Weight (kg) :	0,900	Weight packaging (kg) :	1,000
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



25/08/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference HDPHON13W