

## DORO PHONEEASY110 BLACK VOICE FEEDBACK



Réf. HDPHON27N



## DORO PHONEASY110 BLACK VOICE FEEDBACK



Réf. HDPHON27N

FEATURES			
Medical device	No		
Marquage		Barcode	7322460059538
Guarantee (month):	24	Customs Code	8517110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	19,0
Height (cm) :		Height packaging (cm) :	15,0
Weight (kg) :	0,504	Weight packaging (kg) :	0,504
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:09 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference HDPHON27N