

DUO DORO PHONEEASY110 BLACK VOICE FEEDBACK



Réf. HDPHON29N



DUO DORO PHONEEASY110 BLACK VOICE FEEDBACK



Réf. HDPHON29N

FEATURES			
Medical device	No		
Marquage	Barcode	7322460062194	
Guarantee (month):	24	Customs Code	8517110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	10,0
Length (cm) :	19,0	Length packaging (cm) :	19,0
Height (cm) :	9,8	Height packaging (cm) :	18,0
Weight (kg) :	0,610	Weight packaging (kg) :	0,998
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:09 - This article is available on our website : www.gsh-med.fr with the following reference HDPHON29N